



TGFBR2 (Phospho Tyr284) rabbit pAb

Catalog No	YP-Ab-13126
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB; ELISA
Gene Name	TGFBR2
Protein Name	TGFBR2 (Phospho Tyr284)
Immunogen	Synthesized peptide derived from human TGFBR2 (Phospho Tyr284)
Specificity	This antibody detects endogenous levels of Human,Mouse,Rat TGFBR2 (Phospho Tyr284)
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:1000-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TGF-beta receptor type-2 (TGFR-2;EC 2.7.11.30;TGF-beta type II receptor;Transforming growth factor-beta receptor type II;TGF-beta receptor type II;TbetaR-II)
Observed Band	130
Cell Pathway	Cell membrane ; Single-pass type I membrane protein . Membrane raft .; [Isoform 3]: Secreted .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images