



Cytokeratin 17 (ABT313) Mouse mAb

Catalog No	YP-Ab-15077
Isotype	IgG
Reactivity	Human
Applications	IHC, WB
Gene Name	KRT17
Protein Name	39.1;CK 17;CK-17;Cytokeratin-17;K17;K1C17_HUMAN;Keratin 17;keratin 17 epitope S1;keratin 17 epitope S2;keratin 17 epitope S4;Keratin 17, type I;Keratin;Keratin type I cytoskeletal 17;keratin, type i cytoskeletal 17 [version 1];Keratin-17;KRT17;PC;PC2;PCHC1;type I cytoskeletal 17
Immunogen	Synthesized peptide derived from human Cytokeratin 17
Specificity	The antibody can specifically recognize human CK17 protein. In western blotting of Hela and A431 cell lysates, the antibody can label a 46 kDa band corresponding to CK17.
Formulation	PBS, pH7.2, 0.03% Porcolin 300, containing stabilizing protein
Source	Monoclonal Mouse IgG2a, kappa
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:200-400, WB: Use a concentration of 1 - 5 µg/ml
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	39.1;CK 17;CK-17;Cytokeratin-17;K17;K1C17_HUMAN;Keratin 17;keratin 17 epitope S1;keratin 17 epitope S2;keratin 17 epitope S4;Keratin 17, type I;Keratin;Keratin type I cytoskeletal 17;keratin, type i cytoskeletal 17 [version 1];Keratin-17;KRT17;PC;PC2;PCHC1;type I cytoskeletal 17
Observed Band	
Cell Pathway	Cytoplasmic, Membranous
Tissue Specificity	Tonsil/ Prostate
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images