



Calcitonin (ABT403) Mouse mAb

Catalog No	YP-Ab-15118
Isotype	IgG
Reactivity	Human
Applications	IHC
Gene Name	CALCA CALC1
Protein Name	Alpha CGRP;alpha type CGRP;CALC 1;CALC A;CALC1;calcitonin 1;calcitonin gene related peptide 1;calcitonin gene related peptide I;calcitonin related polypeptide alpha;CGRP 1;CGRP;CGRP I;CGRP1;CT;katacalcin;KC
Immunogen	Synthesized peptide derived from human Calcitonin
Specificity	The antibody can specifically recognize human Calcitonin protein.
Formulation	PBS, pH7.2, 0.03% Porcolin 300, containing stabilizing protein
Source	Monoclonal Mouse IgG1, Kappa
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:200-400,
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Alpha CGRP;alpha type CGRP;CALC 1;CALC A;CALC1;calcitonin 1;calcitonin gene related peptide 1;calcitonin gene related peptide I;calcitonin related polypeptide alpha;CGRP 1;CGRP;CGRP I;CGRP1;CT;katacalcin;KC
Observed Band	
Cell Pathway	Cytoplasmic
Tissue Specificity	Thyroid/ Medullarythyroidcarcinoma
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images