



Cytokeratin 13 (ABT151) Mouse mAb

Catalog No	YP-Ab-15150
Isotype	IgG
Reactivity	Human
Applications	IHC, WB
Gene Name	KRT13
Protein Name	47 kDa cytokeratin;CK-13;CK13;Cytokeratin 13;Cytokeratin-13;K13;K1C13_HUMAN;Ka13;Keratin 13;Keratin;keratin type I cytoskeletal 13;Keratin-13;Krt-1.13;Krt1-13;KRT13;MGC161462;MGC3781;type I cytoskeletal 13;Type I keratin Ka13;WSN2
Immunogen	Synthesized peptide derived from human Cytokeratin 13
Specificity	The antibody can specifically recognize human CK13 protein. In western blotting of A431 cell lysate, the antibody can label a 50 kDa band corresponding to CK13.
Formulation	PBS, pH7.2, 0.03% Porcolin 300, containing stabilizing protein
Source	Monoclonal Mouse IgG1, kappa
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:200-400, WB: Use a concentration of 1 - 5 µg/ml
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	47 kDa cytokeratin;CK-13;CK13;Cytokeratin 13;Cytokeratin-13;K13;K1C13_HUMAN;Ka13;Keratin 13;Keratin;keratin type I cytoskeletal 13;Keratin-13;Krt-1.13;Krt1-13;KRT13;MGC161462;MGC3781;type I cytoskeletal 13;Type I keratin Ka13;WSN2
Observed Band	
Cell Pathway	Cytoplasmic, Membranous
Tissue Specificity	Tonsil
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images