



Fascin (ABT468) Mouse mAb

Catalog No	YP-Ab-15172
Isotype	IgG
Reactivity	Human
Applications	IHC, WB
Gene Name	FSCN1 FAN1 HSN SNL
Protein Name	55 kDa actin bundling protein; 55 kDa actin-bundling protein; Actin bundling protein; actin bundling protein, 55-KD; FAN 1; FAN1; Fascin 1; Fascin actin bundling protein 1; Fascin; Fascin homolog 1 actin bundling protein (Strongylocentrotus purpuratus); Fascin homolog 1; Fascin, sea urchin, homolog of, 1; Fascin1; FLJ38511; FSCN 1; FSCN1; FSCN1_HUMAN; HSN; p55; Singed (Drosophila) like (sea urchin fascin homolog like); Singed drosophila homolog like; Singed like (fascin homolog sea urchin) (Drosophila); Singed like (fascin homolog sea urchin); Singed like protein; Singed, drosophila, homolog of; Singed-like protein; SNL; Strongylocentrotus purpuratus
Immunogen	Synthesized peptide derived from human Fascin
Specificity	The antibody can specifically recognize human Fascin protein. In western blotting of Hela cell lysate, the antibody can label a 55 kDa band corresponding to Fascin.
Formulation	PBS, pH7.2, 0.03% Porcolin 300, containing stabilizing protein
Source	Monoclonal Mouse IgG1, kappa
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:200-400, WB: Use a concentration of 1 - 5 µg/ml
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	55 kDa actin bundling protein; 55 kDa actin-bundling protein; Actin bundling protein; actin bundling protein, 55-KD; FAN 1; FAN1; Fascin 1; Fascin actin bundling protein 1; Fascin; Fascin homolog 1 actin bundling protein (Strongylocentrotus purpuratus); Fascin homolog 1; Fascin, sea urchin, homolog of, 1; Fascin1; FLJ38511; FSCN 1; FSCN1; FSCN1_HUMAN; HSN; p55; Singed (Drosophila) like (sea urchin fascin homolog like); Singed drosophila homolog like; Singed like (fascin homolog sea urchin) (Drosophila); Singed like (fascin homolog sea urchin); Singed like protein; Singed, drosophila, homolog of; Singed-like protein; SNL; Strongylocentrotus purpuratus
Observed Band	
Cell Pathway	Cytoplasmic

**Tissue Specificity**

Tonsil

Function**Background****matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images