



Glycophorin A (ABT334) Mouse mAb

Catalog No	YP-Ab-15178
Isotype	IgG
Reactivity	Human
Applications	IHC
Gene Name	GYPA GPA
Protein Name	AI853584;Blood group--MN locus;CD_antigen=CD235a;CD235a;GLPA_HUMAN;Glycophorin A (MNS blood group);Glycophorin A;Glycophorin A, included;Glycophorin-A;GlycophorinA;GPA;GPErik;GpMiIII;GPSAT;GYPA;GYPA , included;HGpMiIII;HGpMiV;HGpMiX;HGpMiXI;HGpSta(C);MN;MN sialoglycoprotein;MNS;PAS-2;PAS2;Sialoglycoprotein alpha
Immunogen	Synthesized peptide derived from human Glycophorin A
Specificity	The antibody can specifically recognize human Glycophorin A protein.
Formulation	PBS, pH7.2, 0.03% Porcolin 300, containing stabilizing protein
Source	Monoclonal Mouse IgG3, Kappa
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:200-400,
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AI853584;Blood group--MN locus;CD_antigen=CD235a;CD235a;GLPA_HUMAN;Glycophorin A (MNS blood group);Glycophorin A;Glycophorin A, included;Glycophorin-A;GlycophorinA;GPA;GPErik;GpMiIII;GPSAT;GYPA;GYPA , included;HGpMiIII;HGpMiV;HGpMiX;HGpMiXI;HGpSta(C);MN;MN sialoglycoprotein;MNS;PAS-2;PAS2;Sialoglycoprotein alpha
Observed Band	
Cell Pathway	Membranous
Tissue Specificity	Spleen
Function	
Background	

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images