



GCDFP-15 (ABT480) Mouse mAb

Catalog No	YP-Ab-15179
Isotype	IgG
Reactivity	Human
Applications	IHC
Gene Name	PIP GCDFP15 GPIP4
Protein Name	BRST-2;BRST2;GCDFP-15;GCDGP15;gp17;GPIP4;Gross cystic disease fluid protein 15;PIP;PIP_HUMAN;Prolactin induced protein;Prolactin inducible protein precursor;Prolactin-induced protein;Prolactin-inducible protein;SABP;Secretory actin binding protein;Secretory actin-binding protein
Immunogen	Synthesized peptide derived from human GCDFP-15
Specificity	The antibody can specifically recognize human GCDFP-15 protein.
Formulation	PBS, pH7.2, 0.03% Porcolin 300, containing stabilizing protein
Source	Monoclonal Mouse IgG1, Kappa
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:200-400,
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BRST-2;BRST2;GCDFP-15;GCDGP15;gp17;GPIP4;Gross cystic disease fluid protein 15;PIP;PIP_HUMAN;Prolactin induced protein;Prolactin inducible protein precursor;Prolactin-induced protein;Prolactin-inducible protein;SABP;Secretory actin binding protein;Secretory actin-binding protein
Observed Band	
Cell Pathway	Cytoplasmic
Tissue Specificity	Breast
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images