



GST-Pi (ABT346) Mouse mAb

Catalog No	YP-Ab-15181
Isotype	IgG
Reactivity	Human
Applications	IHC, WB
Gene Name	GSTP1 FAEES3 GST3
Protein Name	Deafness; Deafness X-linked 7; DFN7; FAEES3; Fatty Acid Ethyl Ester Synthase III; Glutathione S Transferase 3; Glutathione S Transferase Pi; Glutathione S-transferase P; Glutathione S-transferase pi 1; GST class-pi; GST3; GSTP; Gstp1; GSTP1-1; GSTP1_HUMAN; PI; X linked 7
Immunogen	Synthesized peptide derived from human GST-Pi
Specificity	The antibody can specifically recognize human GST-Pi protein. In western blotting of Hela cell lysate, the antibody can label a 26 kDa band corresponding to GST-Pi.
Formulation	PBS, pH7.2, 0.03% Porcolin 300, containing stabilizing protein
Source	Monoclonal Mouse IgG2a, Kappa
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:200-400, WB: Use a concentration of 1 - 5 µg/ml
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Deafness; Deafness X-linked 7; DFN7; FAEES3; Fatty Acid Ethyl Ester Synthase III; Glutathione S Transferase 3; Glutathione S Transferase Pi; Glutathione S-transferase P; Glutathione S-transferase pi 1; GST class-pi; GST3; GSTP; Gstp1; GSTP1-1; GSTP1_HUMAN; PI; X linked 7
Observed Band	
Cell Pathway	Cytoplasmic
Tissue Specificity	Esophageal squamous cell carcinoma
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images