



MCM2 (ABT195) Mouse mAb

Catalog No	YP-Ab-15183
Isotype	IgG
Reactivity	Human
Applications	IHC, WB
Gene Name	MCM2 BM28 CCNL1 CDCL1 KIAA0030
Protein Name	BM28;CCNL 1;CCNL1;CDC like 1;cdc19;CDCL 1;CDCL1;Cell division cycle like 1;Cyclin like 1;D3S3194;DNA replication licensing factor MCM2;KIAA0030;MCM 2;MCM2;MCM2 minichromosome maintenance deficient 2 mitotin;MCM2 minichromosome maintenance deficient 2 mitotin (S. cerevisiae);MCM2 minichromosome maintenance deficient 2, mitotin;MCM2_HUMAN;MGC10606;Minichromosome maintenance complex component 2;Minichromosome maintenance deficient 2 (mitotin);Minichromosome maintenance deficient 2 mitotin;Minichromosome maintenance protein 2;Minichromosome maintenance protein 2 homolog;Mitotin;Nuclear protein BM28;OTTHUMP00000216047;OTTHUMP00000216050
Immunogen	Synthesized peptide derived from human MCM2
Specificity	The antibody can specifically recognize human MCM2 protein. In western blotting of HEK293 cell lysate, the antibody can label a 130 kDa band.
Formulation	PBS, pH7.2, 0.03% Porcolin 300, containing stabilizing protein
Source	Monoclonal Mouse IgG1, Kappa
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:200-400, WB: Use a concentration of 1 - 5 µg/ml
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BM28;CCNL 1;CCNL1;CDC like 1;cdc19;CDCL 1;CDCL1;Cell division cycle like 1;Cyclin like 1;D3S3194;DNA replication licensing factor MCM2;KIAA0030;MCM 2;MCM2;MCM2 minichromosome maintenance deficient 2 mitotin;MCM2 minichromosome maintenance deficient 2 mitotin (S. cerevisiae);MCM2 minichromosome maintenance deficient 2, mitotin;MCM2_HUMAN;MGC10606;Minichromosome maintenance complex component 2;Minichromosome maintenance deficient 2 (mitotin);Minichromosome maintenance deficient 2 mitotin;Minichromosome maintenance protein 2;Minichromosome maintenance protein 2 homolog;Mitotin;Nuclear protein BM28;OTTHUMP00000216047;OTTHUMP00000216050

**Observed Band****Cell Pathway** Nuclear**Tissue Specificity** Tonsil**Function****Background****matters needing attention** Avoid repeated freezing and thawing!**Usage suggestions** This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.**Products Images**