



p120 Catenin (ABT474) Mouse mAb

Catalog No	YP-Ab-15194
Isotype	IgG
Reactivity	Human
Applications	IHC, WB
Gene Name	CTNND1 KIAA0384
Protein Name	Cadherin associated Src substrate; Cadherin-associated Src substrate; CAS; Catenin (cadherin associated protein) delta 1; Catenin delta 1; Catenin delta; Catenin delta-1; CTND1_HUMAN; CTNND1; CTNND; CTNND1; delta 1 Catenin; KIAA0384; p120; P120 CAS; p120 catenin; P120 CTN; p120(cas); p120(ctn); P120CAS; P120CTN
Immunogen	Synthesized peptide derived from human p120 Catenin
Specificity	The antibody can specifically recognize human p120 catenin protein. In western blotting of Hela, A431 and HEK293 cell lysates, the antibody can label a band with molecular weight around 100~110 kDa.
Formulation	PBS, pH7.2, 0.03% Porcolin 300, containing stabilizing protein
Source	Monoclonal Mouse IgG2b, Kappa
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:200-400, WB: Use a concentration of 1 - 5 µg/ml
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Cadherin associated Src substrate; Cadherin-associated Src substrate; CAS; Catenin (cadherin associated protein) delta 1; Catenin delta 1; Catenin delta; Catenin delta-1; CTND1_HUMAN; CTNND1; CTNND; CTNND1; delta 1 Catenin; KIAA0384; p120; P120 CAS; p120 catenin; P120 CTN; p120(cas); p120(ctn); P120CAS; P120CTN
Observed Band	
Cell Pathway	Cytoplasmic, Membranous
Tissue Specificity	Breast/ Tonsil
Function	
Background	

**matters needing attention**

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images