



MUC2 (ABT-MUC2) mouse mAb

Catalog No	YP-Ab-15362
Isotype	IgG
Reactivity	Human
Applications	IHC-p;IF(paraffin section)
Gene Name	MUC2 SMUC
Protein Name	Mucin-2 (MUC-2) (Intestinal mucin-2)
Immunogen	Synthesized peptide derived from human MUC2
Specificity	This antibody detects endogenous levels of human MUC2. Heat-induced epitope retrieval (HIER) TRIS-EDTA of pH8.0 was highly recommended as antigen repair method in paraffin section
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse/IgG2b, Kappa
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:100-500
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Secreted. In the intestine, secreted into the inner and outer mucus layers. .
Tissue Specificity	Colon, small intestine, colonic tumors, bronchus, cervix and gall bladder.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

