



AMACR (PT1743) mouse mAb

Catalog No	YP-Ab-15374
Isotype	IgG
Reactivity	Human
Applications	IHC, WB
Gene Name	AMACR
Protein Name	Alpha-methylacyl-CoA racemase (EC 5.1.99.4) (2-methylacyl-CoA racemase)
Immunogen	Synthesized peptide derived from human AMACR
Specificity	This antibody detects endogenous levels of human AMACR. Heat-induced epitope retrieval (HIER) TRIS-EDTA of pH8.0 was highly recommended as antigen repair method in paraffin section
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse/IgG1, Kappa
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:100-500, WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	Peroxisome . Mitochondrion .
Tissue Specificity	Aorta, Brain, Cerebellum, Kidney, Liver, PCR rescued clones, Prostate cancer, Sali
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

