



Thyroid transcription factor-1 (TTF1) mouse mAb(ABT-TTF1)

Catalog No	YP-Ab-15458
Isotype	IgG
Reactivity	Human
Applications	IHC-p;IF(paraffin section)
Gene Name	NKX2-1 NKX2A TTF1 TTF1
Protein Name	Thyroid transcription factor-1 (TTF1)
Immunogen	Synthesized peptide derived from human Thyroid transcription factor-1 (TTF1)
Specificity	This antibody detects endogenous levels of human Thyroid transcription factor-1 (TTF1). Heat-induced epitope retrieval (HIER) TRIS-EDTA of pH9.0 was highly recommended as antigen repair method in paraffin section
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse/IgG1, Kappa
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:100-500
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Homeobox protein Nkx-2.1 (Homeobox protein NK-2 homolog A;Thyroid nuclear factor 1;Thyroid transcription factor 1;TTF-1)
Observed Band	
Cell Pathway	Nucleus .
Tissue Specificity	Thyroid and lung.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images