



CK6 mouse mAb(ABT508)

Catalog No	YP-Ab-15583
Isotype	IgG
Reactivity	Human
Applications	IHC, WB
Gene Name	KRT6A K6A KRT6D
Protein Name	CK6
Immunogen	Synthesized peptide derived from human CK6
Specificity	The antibody can specifically recognize human CK6 protein, and shows no cross reaction with CK5. In western blotting of Hela, A431, A549 and HEK293 cell lysates, the antibody can label a 56KDa band corresponding to CK6 only in A431 cell lysates.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.33% sodium azide.
Source	Monoclonal, Mouse/IgG1, kappa
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:100-500, IF 1:100-500
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Keratin, type II cytoskeletal 6A (Cytokeratin-6A; CK-6A; Cytokeratin-6D; CK-6D; Keratin-6A; K6A; Type-II keratin Kb6; allergen Hom s 5)
Observed Band	
Cell Pathway	Tonsil, Squamous cell lung carcinoma
Tissue Specificity	Expressed in the corneal epithelium (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images