



Gross cystic disease fluid protein-15 (GCDFP-15) mouse mAb(ABT469)

Catalog No	YP-Ab-15593
Isotype	IgG
Reactivity	Human
Applications	IHC-p;IF(paraffin section)
Gene Name	PIP GCDFP15 GPII4
Protein Name	Gross cystic disease fluid protein-15 (GCDFP-15)
Immunogen	Synthesized peptide derived from human Gross cystic disease fluid protein-15 (GCDFP-15)
Specificity	The antibody can specifically recognize human GCDFP-15 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.88% sodium azide.
Source	Monoclonal, Mouse/IgG1, kappa
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:100-500, IF 1:100-500
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Prolactin-inducible protein (Gross cystic disease fluid protein 15;GCDFP-15;Prolactin-induced protein;Secretory actin-binding protein;SABP;gp17)
Observed Band	
Cell Pathway	Secreted.
Tissue Specificity	Expressed in pathological conditions of the mammary gland and in several exocrine tissues, such as the lacrimal, salivary, and sweat glands.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images