



TCL1 mouse mAb(ABT431)

Catalog No	YP-Ab-15625
Isotype	IgG
Reactivity	Human
Applications	IHC-p;IF(paraffin section)
Gene Name	TCL1A TCL1
Protein Name	TCL1
Immunogen	Synthesized peptide derived from human TCL1
Specificity	The antibody can specifically recognize human TCL1 protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.152% sodium azide.
Source	Monoclonal, Mouse/IgG1, Kappa
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:100-500, IF 1:100-500
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	T-cell leukemia/lymphoma protein 1A (Oncogene TCL-1;Oncogene TCL1;Protein p14 TCL1)
Observed Band	
Cell Pathway	Cytoplasm . Nucleus . Microsome . Endoplasmic reticulum . Microsomal fraction.
Tissue Specificity	Restricted in the T-cell lineage to immature thymocytes and activated peripheral lymphocytes. Preferentially expressed early in T- and B-lymphocyte differentiation.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images