



Actin, muscle specific(MSA) mouse mAb(ABT024)

Catalog No	YP-Ab-15632
Isotype	IgG
Reactivity	Human
Applications	IHC-p;IF(paraffin section)
Gene Name	
Protein Name	Actin, muscle specific(MSA)
Immunogen	Synthesized peptide derived from human Actin, muscle specific(MSA)
Specificity	The antibody can specifically recognize proteins from smooth muscle, myocardium and skeletal muscle α actin and smooth muscle derived γ actin, an β actin or non smooth muscle derived γ actin does not respond.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.160% sodium azide.
Source	Monoclonal, Mouse/IgG2b, Kappa
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:100-500, IF 1:100-500
Concentration	1 mg/ml
Purity	$\geq 90\%$
Storage Stability	-20°C/1 year
Synonyms	ERROR 42
Observed Band	
Cell Pathway	Cytoplasm, cytoskeleton.
Tissue Specificity	Cytoplasmic
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

