



CD68 (ABT458) Mouse mAb

Catalog No	YP-Ab-15649
Isotype	IgG
Reactivity	Human
Applications	IHC
Gene Name	CD68
Protein Name	CD 68;CD68;CD68 antigen;CD68 molecule;CD68_HUMAN;DKFZp686M18236;gp11;Gp110;LAMP4;Macrophage antigen CD68 (microsialin);MACROPHAGE ANTIGEN CD68;Macrosialin;SCARD1;Scavenger receptor class D member 1
Immunogen	Synthesized peptide derived from human CD68
Specificity	The antibody can specifically recognize human CD68 protein.
Formulation	PBS, pH7.2, 0.03% Porcolin 300, containing stabilizing protein
Source	Monoclonal Mouse IgG1, Kappa
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:200-400,
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CD 68;CD68;CD68 antigen;CD68 molecule;CD68_HUMAN;DKFZp686M18236;gp11;Gp110;LAMP4;Macrophage antigen CD68 (microsialin);MACROPHAGE ANTIGEN CD68;Macrosialin;SCARD1;Scavenger receptor class D member 1
Observed Band	
Cell Pathway	Cytoplasmic, Membranous
Tissue Specificity	Tonsil/ Appendix
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images