



CK16 mouse mAb(ABT312)

Catalog No	YP-Ab-15655
Isotype	IgG
Reactivity	Human
Applications	IHC-p;IF(paraffin section)
Gene Name	
Protein Name	CK16
Immunogen	Synthesized peptide derived from human CK16
Specificity	The antibody can specifically recognize human CK16 protein, and shows no cross reaction with CK1, 5, 6, 7, 8,10, 13, 14, 15, 19, 20.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.73% sodium azide.
Source	Monoclonal, Mouse/IgG1, kappa
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:100-500, IF 1:100-500
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Keratin, type I cytoskeletal 16 (Cytokeratin-16;CK-16;Keratin-16;K16)
Observed Band	
Cell Pathway	Tonsil, Cutaneous squamous cell carcinoma
Tissue Specificity	Expressed in the corneal epithelium (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

