



Factor XIIIa mouse mAb(ABT516)

Catalog No	YP-Ab-15662
Isotype	IgG
Reactivity	Human
Applications	IHC, WB
Gene Name	F13A1 F13A
Protein Name	Factor XIIIa
Immunogen	Synthesized peptide derived from human Factor XIIIa
Specificity	The antibody can specifically recognize human Factor XIIIa protein. In western blotting of A431 cell lysate, there shows a 100KDa band detected by the antibody.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.16% sodium azide.
Source	Monoclonal, Mouse/IgG3, Kappa
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:100-500, IF 1:100-500
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Coagulation factor XIII A chain (Coagulation factor XIIIa; EC 2.3.2.13; Protein-glutamine gamma-glutamyltransferase A chain; Transglutaminase A chain)
Observed Band	
Cell Pathway	Cytoplasm. Secreted . Secreted into the blood plasma. Cytoplasmic in most tissues, but also secreted in the blood plasma.
Tissue Specificity	Cytoplasmic
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images