



Glycophorin A mouse mAb(ABT491)

Catalog No	YP-Ab-15666
Isotype	IgG
Reactivity	Human
Applications	IHC-p;IF(paraffin section)
Gene Name	GYPA GPA
Protein Name	Glycophorin A
Immunogen	Synthesized peptide derived from human Glycophorin A
Specificity	The antibody can specifically recognize human Glycophorin A protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.35% sodium azide.
Source	Monoclonal, Mouse/IgG3, kappa
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:100-500, IF 1:100-500
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Glycophorin-A (MN sialoglycoprotein;PAS-2;Sialoglycoprotein alpha;CD antigen CD235a)
Observed Band	
Cell Pathway	Cell membrane ; Single-pass type I membrane protein . Appears to be colocalized with SLC4A1.
Tissue Specificity	Cytoplasmic
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

