



Lysozyme mouse mAb(ABT342)

Catalog No	YP-Ab-15673
Isotype	IgG
Reactivity	Human
Applications	IHC-p;IF(paraffin section)
Gene Name	LYZ LYM
Protein Name	Lysozyme
Immunogen	Synthesized peptide derived from human Lysozyme
Specificity	The antibody can specifically recognize human Lysozyme protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.157% sodium azide.
Source	Monoclonal, Mouse/IgG3, Kappa
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:100-500, IF 1:100-500
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Lysozyme C (EC 3.2.1.17;1,4-beta-N-acetylmuramidase C)
Observed Band	
Cell Pathway	Secreted.
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

