



Vimentin (ABT283) Mouse mAb

Catalog No	YP-Ab-15685
Isotype	IgG
Reactivity	Human
Applications	IHC, WB
Gene Name	VIM
Protein Name	CTRCT30; Epididymis luminal protein 113; FLJ36605; HEL113; VIM; VIME_HUMAN; Vimentin
Immunogen	Synthesized peptide derived from human Vimentin
Specificity	The antibody can specifically recognize human Vimentin protein. In western blotting of Hela cell lysate, the antibody can label a 54 kDa band corresponding to Vimentin.
Formulation	PBS, pH7.2, 0.03% Porcolin 300, containing stabilizing protein
Source	Monoclonal Mouse IgG1, kappa
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:200-400, WB: Use a concentration of 1 - 5 µg/ml
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CTRCT30; Epididymis luminal protein 113; FLJ36605; HEL113; VIM; VIME_HUMAN; Vimentin
Observed Band	
Cell Pathway	Cytoplasmic
Tissue Specificity	Tonsil/ Liver
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images