



Myoglobin (ABT352) Mouse mAb

Catalog No	YP-Ab-15711
Isotype	IgG
Reactivity	Human
Applications	IHC
Gene Name	MB
Protein Name	MB;MGC13548;MYG_HUMAN;Myoglobin
Immunogen	Synthesized peptide derived from human Myoglobin
Specificity	The antibody can specifically recognize human Myoglobin protein.
Formulation	PBS, pH7.2, 0.03% Porcolin 300, containing stabilizing protein
Source	Monoclonal Mouse IgG2b, Kappa
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:200-400,
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MB;MGC13548;MYG_HUMAN;Myoglobin
Observed Band	
Cell Pathway	Cytoplasmic
Tissue Specificity	Skeletal muscle
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

