



EGF Monoclonal Antibody

Catalog No	YP-Ab-15830
Isotype	IgG
Reactivity	Human
Applications	WB;IF/ICC;ELISA
Gene Name	EGF
Protein Name	Pro-epidermal growth factor
Immunogen	Purified recombinant fragment of human EGF expressed in E. Coli.
Specificity	EGF Monoclonal Antibody detects endogenous levels of EGF protein.
Formulation	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Source	Monoclonal, Mouse
Purification	Affinity purification
Dilution	Western Blot: 1/500 - 1/2000. Immunofluorescence: 1/200 - 1/1000. ELISA: 1/10000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EGF; Pro-epidermal growth factor; EGF
Observed Band	
Cell Pathway	Membrane; Single-pass type I membrane protein.
Tissue Specificity	Expressed in kidney, salivary gland, cerebrum and prostate.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

