



Human EGF mouse mAb

Catalog No	YP-Ab-15858
Isotype	IgG
Reactivity	Human
Applications	ELISA
Gene Name	egf
Protein Name	
Immunogen	Purified recombinant Human EGF protein fragments expressed in E.coli.
Specificity	This antibody detects Human EGF.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Monoclonal, Mouse
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Dilution	ELISA 1:10000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AI790464 ;Beta urogastrone ;Epidermal Growth Factor ;HOMG4 ;OTTHUMP00000219721 ;OTTHUMP00000219722 ;Pro epidermal growth factor ;Pro epidermal growth factor precursor (EGF) ;URG ;Urogastrone .
Observed Band	
Cell Pathway	Membrane; Single-pass type I membrane protein.
Tissue Specificity	Expressed in kidney, salivary gland, cerebrum and prostate.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images