



Anti-Hemoglobin Beta F Monoclonal Mouse Antibody for Lateral Flow Assay Application

Catalog No YP-Ab-17129

Isotype IgG

Reactivity

Applications

Gene Name

Protein Name

Immunogen

Specificity

Formulation

Source

Purification

Dilution

Concentration 1 mg/ml

Purity $\geq 90\%$

Storage Stability $-20^{\circ}\text{C}/1$ year

Synonyms

Observed Band

Cell Pathway

Tissue Specificity

Function

Background

matters needing attention Avoid repeated freezing and thawing!

Usage suggestions This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

