



Recombinant Anti-CD27 Antibody, Rabbit Monoclonal

Catalog No	YP-Ab-17130
Isotype	IgG
Reactivity	
Applications	IHC-P(FAQProtocol)
Gene Name	交叉反应Reacts with: Mouse
Protein Name	
Immunogen	Recombinant Mouse CD27/TNFRSF7 Protein
Specificity	Mouse CD27
Formulation	
Source	
Purification	
Dilution	
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

