



HMGCRC Rabbit mAb

Catalog No	YP-Ab-17137
Isotype	IgG
Reactivity	Human, Mouse
Applications	WB, IP
Gene Name	HMGCRC
Protein Name	LDLCQ3
Immunogen	A synthetic peptide of human HMGCRC
Specificity	Research Field: Cardiovascular
Formulation	Conjugation: Unconjugated
Source	Host: Rabbit
Purification	Clonality: Monoclonal Antibody
Dilution	WB: 1/500-1/1000 IP: 1/20
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Purification: Affinity Purified
Observed Band	97 kDa
Cell Pathway	Modification: Unmodified
Tissue Specificity	
Function	Buffer System: 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Background	Background: Transmembrane glycoprotein that is the rate-limiting enzyme in cholesterol biosynthesis as well as in the biosynthesis of nonsterol isoprenoids that are essential for normal cell function including ubiquinone and geranylgeranyl proteins.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

