



Adipose Triglyceride Lipase Rabbit mAb

Catalog No	YP-Ab-17139
Isotype	IgG
Reactivity	Human, Mouse, Rat
Applications	WB
Gene Name	PNPLA2
Protein Name	Adipose Triglyceride Lipase Rabbit mAb
Immunogen	A synthesized peptide derived from human ATGL / PNPLA2
Specificity	
Formulation	
Source	Monoclonal Antibody
Purification	Affinity Purified Conjugation Unconjugated Modification Unmodified Form Liquid
Dilution	WB: 1/500-1/1000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Desnutrin; plpl; plpl2; Pnpla2; TTS 2.2; TTS2; TTS2.2; ZETA
Observed Band	
Cell Pathway	
Tissue Specificity	Buffer System (缓冲系统): Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Function	Host: Rabbit
Background	Gene ID: 57104
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

