



N-cadherin Monoclonal Antibody

Catalog No	YP-Ab-17672
Isotype	IgG
Reactivity	CDH2; CDHN; NCAD; Cadherin-2; CDw325; Neural cadherin; N-cadherin; CD antigen CD325
Applications	
Gene Name	WB;IHC;IF;FCM;ELISA
Protein Name	>>Cell adhesion molecules;>>Arrhythmogenic right ventricular cardiomyopathy
Immunogen	CDH2
Specificity	Cadherin-2
Formulation	Purified recombinant fragment of human N-cadherin expressed in E. Coli.
Source	N-cadherin Monoclonal Antibody detects endogenous levels of N-cadherin protein.
Purification	Monoclonal, Mouse
Dilution	Ascitic fluid containing 0.03% sodium azide,0.5% BSA, 50%glycerol.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

