



GGT5 (heavy chain, Cleaved-Gly387) rabbit pAb

| | |
|----------------------------------|--|
| Catalog No | YP-Ab-04383 |
| Isotype | IgG |
| Reactivity | Human;Mouse;Rat |
| Applications | WB; ELISA |
| Gene Name | GGT5 GGTLA1 |
| Protein Name | GGT5 (heavy chain, Cleaved-Gly387) |
| Immunogen | Synthesized peptide derived from human GGT5 (heavy chain, Cleaved-Gly387) |
| Specificity | This antibody detects endogenous levels of Human,Mouse,Rat GGT5 (heavy chain, Cleaved-Gly387, protein was cleaved amino acid sequence between 387-388) |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Source | Polyclonal, Rabbit,IgG |
| Purification | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |
| Dilution | WB 1:1000-2000 ELISA 1:5000-20000 |
| Concentration | 1 mg/ml |
| Purity | ≥90% |
| Storage Stability | -20°C/1 year |
| Synonyms | Gamma-glutamyltransferase 5 (GGT 5;EC 2.3.2.2;Gamma-glutamyl transpeptidase-related enzyme;GGT-rel;Gamma-glutamyltransferase-like activity 1;Gamma-glutamyltranspeptidase 5;Glutathione hydrolase 5;EC 3.4.19.13;Leukotriene-C4 hydrolase;EC 3.4.19.14) [Cleaved into: Gamma-glutamyltransferase 5 heavy chain; Gamma-glutamyltransferase 5 light chain] |
| Observed Band | 42 65kD |
| Cell Pathway | Membrane ; Single-pass type II membrane protein . |
| Tissue Specificity | |
| Function | |
| Background | |
| matters needing attention | Avoid repeated freezing and thawing! |

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images