




iHog Monoclonal Antibody

Catalog No	YP-Ab-04447
Isotype	IgG
Reactivity	Drosophilamelanogaster
Applications	WB;ELISA
Gene Name	
Protein Name	
Immunogen	Purified recombinant fragment of human iHog expressed in E. Coli.
Specificity	iHog Monoclonal Antibody detects iHog protein.
Formulation	Ascitic fluid containing 0.03% sodium azide, 0.5% BSA, 50% glycerol.
Source	Monoclonal, Mouse
Purification	Affinity purification
Dilution	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western Blot analysis using iHog Monoclonal Antibody against IHOG-hlgGFc transfected HEK293 cell lysate.