



# His-Tag Monoclonal Antibody(4E6), FITC Conjugated

|                                  |                                                                                                                                             |
|----------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Catalog No</b>                | YP-Ab-04594                                                                                                                                 |
| <b>Isotype</b>                   | IgG                                                                                                                                         |
| <b>Reactivity</b>                | Species independent                                                                                                                         |
| <b>Applications</b>              | WB;IP;IF                                                                                                                                    |
| <b>Gene Name</b>                 | His-Tag                                                                                                                                     |
| <b>Protein Name</b>              | His Tag                                                                                                                                     |
| <b>Immunogen</b>                 |                                                                                                                                             |
| <b>Specificity</b>               | His-Tag Monoclonal Antibody(4E6) FITC conjugated specially designed for your WB or IHC analysis.                                            |
| <b>Formulation</b>               | Liquid in PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol.                                                      |
| <b>Source</b>                    | Monoclonal, Mouse IgG                                                                                                                       |
| <b>Purification</b>              | The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.                                  |
| <b>Dilution</b>                  | Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: IF: 1:1000. |
| <b>Concentration</b>             | 1 mg/ml                                                                                                                                     |
| <b>Purity</b>                    | ≥90%                                                                                                                                        |
| <b>Storage Stability</b>         | -20°C/1 year                                                                                                                                |
| <b>Synonyms</b>                  | His                                                                                                                                         |
| <b>Observed Band</b>             |                                                                                                                                             |
| <b>Cell Pathway</b>              |                                                                                                                                             |
| <b>Tissue Specificity</b>        |                                                                                                                                             |
| <b>Function</b>                  |                                                                                                                                             |
| <b>Background</b>                |                                                                                                                                             |
| <b>matters needing attention</b> | Avoid repeated freezing and thawing!                                                                                                        |
| <b>Usage suggestions</b>         | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.           |



## Products Images

