



# His-Tag Monoclonal Antibody(4E6)

|                           |   |
|---------------------------|---|
| Catalog No                | YP-Ab-04742   |
| Isotype                   | IgG   |
| Reactivity                | Species independent   |
| Applications              | WB  |
| Gene Name                 | His-Tag   |
| Protein Name              | His Tag   |
| Immunogen                 | Synthetic peptide:HHHHHH conjugated to KLH  |
| Specificity               | The antibody detects C-terminal, internal, and N-terminal His-tag fusion proteins.  |
| Formulation               | PBS, pH 7.4, containing 0.5%BSA, 0.02% sodium azide as Preservative and 50% Glycerol.   |
| Source                    | Monoclonal, Mouse   |
| Purification              | The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.                        |
| Dilution                  | WB: 1:3000-10000  |
| Purity                    | ≥90%  |
| Storage Stability         | -20°C/1 year  |
| matters needing attention | Avoid repeated freezing and thawing!  |
| Usage suggestions         | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel. |

## Products Images