



c-Jun Mouse mAb

Catalog No	YP-mAb-17706
Isotype	IgG
Reactivity	Human,Mouse,Monkey
Applications	WB
Gene Name	c-Jun
Protein Name	
Immunogen	Purified recombinant fragment of human c-Jun expressed in E. Coli.
Specificity	
Formulation	
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AP1; AP-1; c-Jun; Jun
Observed Band	
Cell Pathway	
Tissue Specificity	
Function	
Background	This gene is the putative transforming gene of avian sarcoma virus 17. It encodes a protein which is highly similar to the viral protein, and which interacts directly with specific target DNA sequences to regulate gene expression. This gene is intronless and is mapped to 1p32-p31, a chromosomal region involved in both translocations and deletions in human malignancies.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images