



RPUSD2 Rabbit pAb

| | |
|---------------------------|---|
| Catalog No | YP-Ab-18647 |
| Isotype | IgG |
| Reactivity | Human, Mouse |
| Applications | WB |
| Gene Name | RPUSD2 C15orf19 |
| Protein Name | RNA pseudouridylate synthase domain-containing protein 2 |
| Immunogen | Synthesized peptide derived from human RPUSD2 |
| Specificity | This antibody detects endogenous levels of RPUSD2 at Human, Mouse |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Source | |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Dilution | WB 1:500-2000 |
| Concentration | 1 mg/ml |
| Purity | ≥90% |
| Storage Stability | -20°C/1 year |
| Synonyms | |
| Observed Band | 60kD |
| Cell Pathway | |
| Tissue Specificity | |
| Function | Pseudouridine synthase that catalyzes pseudouridylation of mRNAs. |
| Background | |
| matters needing attention | Avoid repeated freezing and thawing! |
| Usage suggestions | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel. |



Products Images

