



HEYL Rabbit pAb

Catalog No	YP-Ab-18658
Isotype	IgG
Reactivity	Human, Mouse
Applications	WB
Gene Name	HEYL BHLHB33 HRT3
Protein Name	Hairy/enhancer-of-split related with YRPW motif-like protein (hHeyL) (Class B basic helix-loop-helix protein 33) (bHLHb33) (Hairy-related transcription factor 3) (HRT-3) (hHRT3)
Immunogen	Synthesized peptide derived from human HEYL
Specificity	This antibody detects endogenous levels of HEYL at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	36kD
Cell Pathway	Nucleus .
Tissue Specificity	
Function	Downstream effector of Notch signaling which may be required for cardiovascular development (By similarity). Transcriptional repressor which binds preferentially to the canonical E box sequence 5'-CACGTG-3' (By similarity). Represses transcription by the cardiac transcriptional activators GATA4 and GATA6.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images