



BEND5 Mouse mAb

Catalog No	YP-mAb-18690
Isotype	IgG
Reactivity	Human, Mouse
Applications	WB
Gene Name	BEND5 C1orf165
Protein Name	BEN domain-containing protein 5
Immunogen	Synthesized peptide derived from human BEND5
Specificity	This antibody detects endogenous levels of BEND5 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from mouse antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	46kD
Cell Pathway	
Tissue Specificity	
Function	Acts as a transcriptional repressor .
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

