



C2orf79/PTRHD1 Rabbit pA

Catalog No	YP-Ab-18709
Isotype	IgG
Reactivity	Human, Mouse
Applications	WB
Gene Name	PTRHD1 C2orf79
Protein Name	Putative peptidyl-tRNA hydrolase PTRHD1 (Peptidyl-tRNA hydrolase domain-containing protein 1)
Immunogen	Synthesized peptide derived from human C2orf79/PTRHD1
Specificity	This antibody detects endogenous levels of C2orf79/PTRHD1 at Human, Mouse
Formulation	
Source	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	15kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

