



# CDX1 Rabbit pAb

<b>Catalog No</b>	YP-Ab-18727
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human,Mouse,Rat
<b>Applications</b>	WB
<b>Gene Name</b>	CDX1
<b>Protein Name</b>	Homeobox protein CDX-1 (Caudal-type homeobox protein 1)
<b>Immunogen</b>	Synthesized peptide derived from human CDX1
<b>Specificity</b>	This antibody detects endogenous levels of CDX1 at Human, Mouse,Rat
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Source</b>	
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Dilution</b>	WB 1:500-2000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	
<b>Observed Band</b>	29kD
<b>Cell Pathway</b>	Nucleus.
<b>Tissue Specificity</b>	Intestinal epithelium.
<b>Function</b>	Plays a role in transcriptional regulation . Involved in activated KRAS-mediated transcriptional activation of PRKD1 in colorectal cancer (CRC) cells . Binds to the PRKD1 promoter in colorectal cancer (CRC) cells . Could play a role in the terminal differentiation of the intestine. Binds preferentially to methylated DNA .
<b>Background</b>	
<b>matters needing attention</b>	Avoid repeated freezing and thawing!
<b>Usage suggestions</b>	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



## Products Images