



# ARIH2 Rabbit pAb

<b>Catalog No</b>	YP-Ab-18737
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human, Mouse
<b>Applications</b>	WB
<b>Gene Name</b>	ARIH2 ARI2 TRIAD1 HT005
<b>Protein Name</b>	E3 ubiquitin-protein ligase ARIH2 (ARI-2) (Protein ariadne-2 homolog) (Triad1 protein)
<b>Immunogen</b>	Synthesized peptide derived from human ARIH2
<b>Specificity</b>	This antibody detects endogenous levels of ARIH2 at Human, Mouse
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Source</b>	
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Dilution</b>	WB 1:500-2000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	
<b>Observed Band</b>	54kD
<b>Cell Pathway</b>	Nucleus . Cytoplasm .
<b>Tissue Specificity</b>	Widely expressed with higher expression in granulocytes.
<b>Function</b>	E3 ubiquitin-protein ligase, which catalyzes ubiquitination of target proteins together with ubiquitin-conjugating enzyme E2 UBE2L3 . Acts as an atypical E3 ubiquitin-protein ligase by working together with cullin-5-RING ubiquitin ligase complex (ECS complex, also named CRL5 complex) and initiating ubiquitination of ECS substrates: associates with ECS complex and specifically mediates addition of the first ubiquitin on ECS targets (By similarity). The initial ubiquitin is then elongated (By similarity). E3 ubiquitin-protein ligase activity is activated upon binding to neddylated form of the ECS complex . Mediates 'Lys-6', 'Lys-48'- and 'Lys-63'-linked polyubiquitination . May play a role in myelopoiesis .
<b>Background</b>	
<b>matters needing attention</b>	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

**Products Images**