



# CCDC40 Mouse pAb

<b>Catalog No</b>	YP-Ab-18754
<b>Isotype</b>	IgG
<b>Reactivity</b>	Human,Mouse
<b>Applications</b>	WB
<b>Gene Name</b>	CCDC40 KIAA1640
<b>Protein Name</b>	Coiled-coil domain-containing protein 40
<b>Immunogen</b>	Synthesized peptide derived from human CCDC40
<b>Specificity</b>	This antibody detects endogenous levels of CCDC40 at Human, Mouse
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Source</b>	
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Dilution</b>	WB 1:500-2000
<b>Concentration</b>	1 mg/ml
<b>Purity</b>	≥90%
<b>Storage Stability</b>	-20°C/1 year
<b>Synonyms</b>	
<b>Observed Band</b>	126kD
<b>Cell Pathway</b>	Cytoplasm . Cell projection, cilium . Localizes to cytoplasm and motile cilium. .
<b>Tissue Specificity</b>	
<b>Function</b>	Required for assembly of dynein regulatory complex (DRC) and inner dynein arm (IDA) complexes, which are responsible for ciliary beat regulation, thereby playing a central role in motility in cilia and flagella . Probably acts together with CCDC39 to form a molecular ruler that determines the 96 nanometer (nm) repeat length and arrangements of components in cilia and flagella (By similarity). Not required for outer dynein arm complexes assembly. Required for axonemal recruitment of CCDC39 .
<b>Background</b>	
<b>matters needing attention</b>	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

**Products Images**