



TNS4 Rabbit pAb

Catalog No	YP-Ab-18772
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	TNS4 CTEN PP14434
Protein Name	Tensin-4 (C-terminal tensin-like protein)
Immunogen	Synthesized peptide derived from human TNS4
Specificity	This antibody detects endogenous levels of TNS4 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	79kD
Cell Pathway	Prostate and placenta. Down regulated in prostate cancer.
Tissue Specificity	
Function	Cell junction, focal adhesion . Cytoplasm, cytoskeleton .
Background	May be involved in cell migration, cartilage development and in linking signal transduction pathways to the cytoskeleton (By similarity). May promote apoptosis, via its cleavage by caspase-3.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

