



ZNF395 Rabbit pAb

Catalog No	YP-Ab-18776
Isotype	IgG
Reactivity	Human
Applications	WB
Gene Name	ZNF395 HDBP2 PBF
Protein Name	Zinc finger protein 395 (HD-regulating factor 2) (HDRF-2) (Huntington disease gene regulatory region-binding protein 2) (HD gene regulatory region-binding protein 2) (HDBP-2) (Papillomavirus regulatory factor 1) (PRF-1) (Papillomavirus-binding factor)
Immunogen	Synthesized peptide derived from human ZNF395
Specificity	This antibody detects endogenous levels of ZNF395 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	56kD
Cell Pathway	Cytoplasm. Nucleus. May shuttle between nucleus and cytoplasm.
Tissue Specificity	Widely expressed.
Function	Plays a role in papillomavirus genes transcription.
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images