



MTMR4Rabbit pAb

Catalog No	YP-Ab-18780
Isotype	IgG
Reactivity	Human,Mouse
Applications	WB
Gene Name	MTMR4 KIAA0647 ZFYVE11
Protein Name	Myotubularin-related protein 4 (FYVE domain-containing dual specificity protein phosphatase 2) (FYVE-DSP2) (Zinc finger FYVE domain-containing protein 11)
Immunogen	Synthesized peptide derived from human MTMR4
Specificity	This antibody detects endogenous levels of MTMR4 at Human, Mouse
Formulation	
Source	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	131kD
Cell Pathway	Cytoplasm . Membrane ; Peripheral membrane protein . Localized to perinuclear region.
Tissue Specificity	Expressed in brain, heart, kidney, spleen, liver, colon, testis, muscle, placenta, thyroid gland, pancreas, ovary, prostate, skin, peripheral blood, and bone marrow.
Function	Dephosphorylates proteins phosphorylated on Ser, Thr, and Tyr residues and low molecular weight phosphatase substrate para-nitrophenylphosphate. Phosphorylates phosphatidylinositol 3,4,5-trisphosphate (PIP3).
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images