



MOXD1 Rabbit pAb

| | |
|----------------------------------|---|
| Catalog No | YP-Ab-18786 |
| Isotype | IgG |
| Reactivity | Human,Mouse |
| Applications | WB |
| Gene Name | MOXD1 MOX UNQ2493/PRO5780 |
| Protein Name | DBH-like monooxygenase protein 1 (Monooxygenase X) |
| Immunogen | Synthesized peptide derived from human MOXD1 |
| Specificity | This antibody detects endogenous levels of MOXD1 at Human, Mouse |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Source | |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Dilution | WB 1:500-2000 |
| Concentration | 1 mg/ml |
| Purity | ≥90% |
| Storage Stability | -20°C/1 year |
| Synonyms | |
| Observed Band | 67kD |
| Cell Pathway | Endoplasmic reticulum membrane ; Single-pass type I membrane protein . |
| Tissue Specificity | Highly expressed in lung, kidney, brain and spinal cord. |
| Function | |
| Background | |
| matters needing attention | Avoid repeated freezing and thawing! |
| Usage suggestions | This product can be used in immunological reaction related experiments. For more information, please consult technical personnel. |



Products Images