



TRIM15 Rabbit pAb

Catalog No	YP-Ab-18792
Isotype	IgG
Reactivity	Human
Applications	WB
Gene Name	TRIM15 RNF93 ZNF178 ZNFB7
Protein Name	Tripartite motif-containing protein 15 (RING finger protein 93) (Zinc finger protein 178) (Zinc finger protein B7)
Immunogen	Synthesized peptide derived from human TRIM15
Specificity	This antibody detects endogenous levels of TRIM15 at Human
Formulation	
Source	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	51kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

