



sesn2 rabbit pAb

Catalog No	YP-Ab-12582
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	SESN2 SEST2
Protein Name	sesn2
Immunogen	Synthesized peptide derived from human sesn2
Specificity	This antibody detects endogenous levels of sesn2 at Human, Mouse,Rat
Formulation	Liquid in PBS containing 50% glycerol, and 0.74% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Sestrin-2 (Hi95)
Calculated Molecular Weight	53kD
Cell Pathway	Cytoplasm .
Tissue Specificity	Widely expressed.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

