



SMTNL2 Rabbit pAb

Catalog No	YP-Ab-18799
Isotype	IgG
Reactivity	Human,Mouse,Rat
Applications	WB
Gene Name	SMTNL2
Protein Name	Smoothelin-like protein 2
Immunogen	Synthesized peptide derived from human SMTNL2
Specificity	This antibody detects endogenous levels of SMTNL2 at Human, Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	
Observed Band	
Calculated Molecular Weight	51kD
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images